Medion Diagnostics

Search-Cyte® Search-Cyte® Plus Search-Cyte® TCS Reagent Red Blood Cells 3±1%

U.S. License No. 1740

For use in the detection of unexpected antibodies

For in vitro diagnostic use

Summary and Principle

Search-Cyte®, Search-Cyte® Plus and Search-Cyte® TCS are suspensions of group O Reagent Red Blood Cells differing in antigenic configuration for use in antibody screening, and are selected to enable detection of most clinically significant blood group antibodies.

The antigen typings of each donor are provided on the antigenic constitution matrix that accompanies each product. Search-Cyte® Reagent Red Blood Cells may be utilized in any of the commonly accepted antibody detection procedures.^{1,2}

Antibodies react with Red Blood Cells possessing the corresponding antigenic determinants. Antibodies may agglutinate Red Blood Cells in saline, low ionic strength solution (LISS), high protein media or antiglobulin testing.

Reagent

Search-Cyte® Reagent Red Blood Cells I and II are individual suspensions of group O Red Blood Cells from two donors, one rr (cde/cde) and one R₁R₂ (CDe/CDe).

Search-Cyte® Plus Reagent Red Blood Cells I and II are individual suspensions of group O Red Blood Cells from two donors, one R₁R₁ (CDe/CDe) or R₁R^W (CDe/C^WDe) and one R₂R₂ (cDe/cDe).

Search-Cyte® TCS Reagent Red Blood Cells I, II and III are individual suspensions of group O Red Blood Cells from three donors, one R₁R₁ (CDe/CDe) or R₁R^W (CDe/C^WDe), one R₂R₂ (cDe/cDe) and one rr (cde/cde).

All Red Blood Cells are 3±1% suspensions in isotonic medium with added buffers (bicarbonate and phosphate) and preservatives (0.03% neomycin and 0.05% chloramphenicol). The suspending medium contains no ingredients to inhibit complement mediated hemolysis. Frozen/thawed Red Blood Cells may have been used in this product. No U.S. standard of potency. Meets other FDA requirements.

Store at 2–8°C. **Do not freeze.** Resuspend by gently inverting immediately prior to use. Reagent Red Blood Cells are ready to use; washing Reagent Red Blood Cells before use is optional. Indication of deterioration: notable hemolysis which may be caused by microbial contamination or improper handling, darkening of Reagent Red Blood Cells or spontaneous clumping.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

The pipette of the vial contains natural rubber latex which may cause allergic reactions.

Specimen Collection and Preparation

No special preparation of the patient is required prior to specimen collection. Serum from freshly clotted blood is preferred. For optimum test results, serum should be stored at 2–8°C no longer than 48 hours prior to testing;

however, serum may be frozen at -20 to -80°C and tested at a later time if necessary. Plasma samples may be used, however, use of plasma may result in failure to detect complement dependent antibodies due to its low complement activity.^{1,2}

Procedure

Reagent Provided

Search-Cyte® Reagent Red Blood Cells I and II, or Search-Cyte® Plus Reagent Red Blood Cells I and II, or Search-Cyte® TCS Reagent Red Blood Cells I, II and III

Materials Required but Not Provided

1. Disposable test tubes (12 x 75 mm or 10 x 75 mm)
2. Physiologic saline
3. Anti-human globulin
4. Potentiator
5. Coombs control cells
6. Optical aid³
7. Centrifuge (calibrated for 1000 rcf* or 150 rcf*)
8. Water bath or heating block at 37°C
9. Timer
10. Pipets

Procedure Outline

If potentiator is to be used, follow manufacturer's instructions for preparation and testing of samples instead of the procedure suggested below.

1. Place 2 or more drops of the test serum in each properly labeled tube.
2. Add 1 drop of appropriate Reagent Red Blood Cells suspension to each tube. Shake all tubes to mix.
3. If immediate spin testing is desired, centrifuge for 15–20 seconds at approximately 1000 rcf* (1 minute at approximately 150 rcf*) or time appropriate to the calibration of each centrifuge.
4. Gently resuspend Red Blood Cells completely and examine immediately for agglutination or hemolysis. Grade and record results.

Procedure A

Indirect Antiglobulin Test

5. Incubate at 37°C for 15–60 minutes.
6. Immediately centrifuge, examine and interpret as in Steps 3 and 4. Record results of 37°C testing.
7. Fill each tube with physiologic saline added in a forceful stream. Centrifuge to pack Red Blood Cells. Carefully decant supernatant. Shake to resuspend Red Blood Cells.
8. Repeat Step 7 twice for a total of 3 washings.
9. Add anti-human globulin to each tube according to manufacturer's instructions and shake to mix.
10. Centrifuge, examine and interpret as in Steps 3 and 4. Grade and record results.
11. Negative reactions obtained with anti-human globulin should be confirmed by adding 1 drop Coombs control cells (see Coombs Control Cells package insert).

* rcf = 0.0001118 x rotation radius (cm) x rpm²

Procedure B

Saline Room Temperature

- Incubate at room temperature (20–25°C) for 15–30 minutes.
- Centrifuge, examine and interpret as in Steps 3 and 4. Grade and record results.

Note: The reactions should be interpreted immediately after centrifugation due to the possibility of dissociation of the antigen-antibody complex.

Quality Control

Known relatively weak positive serums should be tested with each Search-Cyte® suspension each day Reagent Red Blood Cells are in use. This testing may be conveniently performed using Blood Bank Quality Assurance (BBQA II) Testing Reagents, available from Medion Diagnostics AG.

Use of an autocontrol may be helpful in distinguishing autoantibodies and alloantibodies. If the autocontrol is positive, the serum may contain auto-antibody and further testing may be indicated.²

To ensure proper centrifugation, each individual centrifuge should be calibrated for the specific test procedure being performed. Red Blood Cells should be packed firmly, but negative control Red Blood Cells should resuspend easily.³

Results

Agglutination and/or hemolysis in any of the Search-Cyte® tubes at any phase of the test procedure prior to the addition of Coombs control cells indicates the presence of unexpected antibodies. Such antibodies are usually directed against the known antigens present on the screening Red Blood Cells, but may be directed against an antigen not indicated on the antigenic constitution matrix. The lack of both agglutination and hemolysis in the test procedure indicates the absence of antibodies to antigens contained in the reagent. Agglutination or hemolysis in the autocontrol tube indicates further studies are necessary.¹

Limitations of Procedure

A 15-minute incubation at 37°C may not be adequate to detect some weak blood group antibodies if no potentiating medium is added to the test system.

Low frequency antigens may not always be present on Search-Cyte® Reagent Red Blood Cells, and a double dose of antigen may be required to detect very weakly reacting antibodies; therefore, negative reactions with the screening Red Blood Cells do not always indicate the absence of unexpected antibodies. Because of the high incidence of the *Fy^a* gene in the Black population, it cannot be assumed that the phenotypes *Fy* (a+b-) and *Fy* (a-b+) in Black donors represent homozygous expressions of the *Fy^a* or *Fy^b* genes. Use of an autocontrol may be helpful in distinguishing between autoantibodies and alloantibodies.²

As in all serological tests, such factors as contaminated materials, improper incubation time or temperature, improper centrifugation or improper examination for agglutination may give rise to false test results.

False negative results may occur if

- Red Blood Cells are not properly washed or human globulins are present as contaminants in glassware. These residual globulins will neutralize the globulin-reactive antibodies present in the anti-human globulin serum.
- antibody elutes from Red Blood Cells during incubation or washing.
- Red Blood Cells and/or serum are stored improperly and lose reactivity.
- anti-human globulin is inadvertently omitted.
- Red Blood Cells are improperly centrifuged.
- incubation times and/or temperatures are incorrect for proper Red Blood Cells sensitization.
- resuspension technique is too vigorous to preserve agglutination of weakly sensitized Red Blood Cells.

False positive results may occur if

- Red Blood Cells have microbial contamination.
- Red Blood Cells are improperly centrifuged.
- antibodies to antibiotics or to other ingredients in the Red Blood Cells suspending medium or in the potentiators used are present in the test serum.
- incomplete resuspension may counterfeit agglutination.
- in rare cases, the test serum contains an antibody directed at one of the components of the reagent diluent.

Specific Performance Characteristics

Each lot of Search-Cyte® Reagent Red Blood Cells is carefully prepared to permit detection of antibodies to the selected Red Blood Cells antigens when used as outlined in these procedures.

All antigen typings listed on the antigenic constitution matrix are confirmed using two sources of antiserum except for those indicated on the antigenic constitution matrix enclosed with each lot.

Identified low incidence antigens present are indicated on the antigenic constitution matrix. Direct antiglobulin tests are negative on all Red Blood Cells.

As with all Red Blood Cells, the reactivity of the product may decrease during the dating period. The rate at which antigen reactivity is lost is partially dependent upon individual donor characteristics that are neither controlled nor predictable by the manufacturer. However, if properly stored when not in use, the reagent can be expected to perform as described throughout its dating.

Bibliography

- Mollison P.L., Blood transfusion in clinical medicine. 10th ed. Blackwell Scientific Publications; 1997: Chapter 8.
- Technical manual of the American Association of Blood Banks. 14th ed. 2002, Chapter 18 and 20.
- Ibidem: Method 8.5, p.765.

Warranty

This product is warranted to perform as described in its labeling and in the product literature, and Medion Diagnostics AG disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Medion Diagnostics AG be liable for any consequential damages arising out of the aforesaid express warranty.

Search-Cyte® Search-Cyte® Plus Search-Cyte® TCS Hematíes reactivos 3±1%

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