

Medion Diagnostics AG

Search-Cyte® Plus Reagent Red Blood Cells 3±1%

Name / I.D. _____ Date of birth _____ Date _____

Lot No. **CEP-1105** Exp. Date **22 NOV 08** Result _____ Signature _____

Vial No.	Rh Pheno-type	Donor Code	Special Typings	Rhesus								MNS				P	Lewis		Lutheran		Kell				Duffy		Kidd		Sex linked								
				D	C	E	c	e	f	C ^w	V	M	N	S	s	P ₁	Le ^a	Le ^b	Lu ^a	Lu ^b	K	k	Kp ^a	Js ^a	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Xg ^a								
I	R ₁ R ₁	K836CC		+	+	0	0	+	NT	0	0	+	+	0	+	0	0	+	0	+	+	+	0	0	0	+	+	0	0	+	I						
II	R ₂ R ₂	M2198TR		+	0	+	+	0	NT	0	0	+	+	+	+	+	0	0	+	0	+	0	0	+	+	0	+	0	0	II							
Auto																																					

Manufactured by:
Medion Diagnostics AG, Bonnstrasse 9, CH-3186 Düringen/Switzerland

U.S. License No. 1740

Remarks see reverse side