

Medion Diagnostics AG

Search-Cyte® TCS Reagent Red Blood Cells 3±1%

Name / I.D. _____ Date of birth _____ Date _____

Lot No. CET-506 Exp. Date 22 NOV 08 Result _____ Signature _____

| Vial No. | Rh Phenotype | Donor Code | Special Typings | Rhesus | | | | | | | | MNS | | | | P | Lewis | | Lutheran | | Kell | | | | Duffy | | Kidd | | Sex linked | | | | | | | | |
|----------|-------------------------------|------------|-----------------|--------|---|---|---|---|----|----------------|---|-----|---|---|---|-----------------|-----------------|-----------------|-----------------|-----------------|------|---|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|------------|-----|---|--|--|--|--|-----------------|--|
| | | | | D | C | E | c | e | f | C ^w | V | M | N | S | s | P ₁ | Le ^a | Le ^b | Lu ^a | Lu ^b | K | k | Kp ^a | Js ^a | Fy ^a | Fy ^b | Jk ^a | Jk ^b | | | | | | | | Xg ^a | |
| I | R ₁ R ₁ | M2230CC | | + | + | 0 | 0 | + | NT | 0 | 0 | + | 0 | + | + | + | 0 | + | 0 | + | 0 | + | 0 | 0 | + | 0 | + | + | + | + | I | | | | | | |
| II | R ₂ R ₂ | M1201CC | | + | 0 | + | + | 0 | NT | 0 | 0 | + | 0 | + | 0 | 0 | + | 0 | 0 | + | 0 | + | 0 | 0 | + | + | 0 | + | + | II | | | | | | | |
| III | rr | M1674TR | | 0 | 0 | 0 | + | + | + | 0 | 0 | 0 | + | 0 | + | + ^{vw} | 0 | + | 0 | + | + | + | 0 | 0 | 0 | + | + | 0 | + | III | | | | | | | |

Manufactured by:
Medion Diagnostics AG, Bonnstrasse 9, CH-3186 Dürdingen/Switzerland

U.S. License No. 1740

Remarks see reverse side

Auto